



Customized Solutions

# Yxoss CBR<sup>®</sup>

Customized Bone Regeneration

marketed by

**Geistlich**

the regeneration  
experts



# Yxoss CBR<sup>®</sup> – The future is now

## Discover the first generation of Customized Bone Regeneration (CBR<sup>®</sup>) for complex bone defects

Yxoss CBR<sup>®</sup> has revolutionized bone augmentation by customizing titanium scaffold using a digital workflow to fit the individual anatomy of each patient. Cutting edge CAD/CAM technology is utilized to satisfy patient-specific requirements for complex bone procedures.



This patented process results in a contoured, form-stable, 3-D printed scaffold of the purest titanium based on CT or CBCT images.

Time-consuming impressions, cutting, shaping and adapting are no longer necessary. Sharp edges from cutting conventional meshes are entirely eliminated.

The combination of Yxoss CBR<sup>®</sup> with Geistlich Bio-Gide<sup>®</sup> prevents soft tissue ingrowth and particle migration while adding Geistlich Bio-Oss<sup>®</sup> to the graft provides long-term stability.


Now, with the recent introduction of Yxoss CBR<sup>®</sup> protect, you have even more choice when developing your treatment plan.

It's a solution made just for you.  
**Exactly.**


# Best-in-class biomaterials

Create an unmatched combination

Geistlich Bio-Gide®



Yxoss CBR® classic



OR


Yxoss CBR® protect



Geistlich Bio-Oss®



**EXACTLY**  
like no other.



# Yxoss CBR<sup>®</sup> classic


## A treatment evolution built on quality and precision

Bring a customized approach to your next planned bone augmentation.  
The open structure of Yxoss CBR<sup>®</sup> classic enables periosteal vascularization that is essential for bone regeneration.

Innovative solution for regenerating **complex alveolar bone defects** by taking advantage of CBCT data combined with **3-D printing technology**

For ideal fixation of Yxoss CBR<sup>®</sup>, a fixation system with diameters from 1.3 to 1.5 mm is recommended.

3-D



Optional: Unique **integrated implant positioning** in the surgical planning – **Yxoss CBR<sup>®</sup> Backward**



Opportunity to **reduce surgery time** without complex adaptations



High **stability** and space maintenance



**Easy Removal Design<sup>®</sup>** with pre-defined breaking points for easy removal of Yxoss CBR<sup>®</sup> on re-entry




ReOss<sup>®</sup> calculates the necessary **augmentation volume** for your case planning

# Yxoss CBR<sup>®</sup> protect

## An added innovation designed for easy removal

In certain cases, excessive bone formation in the apical part of the titanium scaffold can hamper its removal due to soft and hard tissue ingrowth. To overcome this challenge Yxoss CBR<sup>®</sup> protect features a dense microporous structure in the apical area which allows an even easier removal.

The periosteal blood supply continues to be promoted in the upper part of the scaffold with the open structure



### Dense microstructure for apical edge zones\*

- › Protects the apical area from soft and hard tissue ingrowth
- › Enables even easier removal




## Biological background

### Periosteal nutrition:

Maturation of the bone is dependent on the vascularization of the periosteal vessels

### No ingrowth of soft and hard tissue:

Higher occlusivity due to denser Yxoss CBR<sup>®</sup> protect structure leads to less intergrowth with tissue cells




\* Manufacturing result may visually differ from rendered microstructure.



Surgery and concept by  
Dr. Philippe Russe  
(Reims, France)

# Horizontal/Vertical Defect (3 teeth gap) – Posterior Mandible



1 CBCT of initial situation showing the mandibular knife edge ridge.

2 Bony situation after full thickness flap.

3 Yxoss CBR®, filled with a mixture (50/50) of autogenous bone and Geistlich Bio-Oss® small granules.

4 Yxoss CBR® placed *in situ*.

5 A Geistlich Bio-Gide® collagen membrane is placed on top, covering the titanium scaffold.

6 Healing at 6 months postoperative: The scaffold is visible under the thin lingual gingiva.

7 Control CBCT at 6 months postoperative.

8 Separation of the Yxoss CBR® into two parts.

9 Implants placed in the bony reconstruction (BEGO Semados® SC).

10 Additional GBR with Geistlich Bio-Oss® alone, Geistlich Bio-Gide®, sutured with horizontal mattress (3/0) and crestal X sutures (6/0).


11 Attached gingiva created by a partial thickness flap around the healing abutments post uncovering of implants after a healing period of 6 months.

12 Provisional bridge at time of ceramic try-in.

# Horizontal/Vertical Defect (3 teeth gap) – Posterior Maxilla



Surgery and concept by  
Dr. Marcus Seiler MSc MSc  
(Filderstadt, Germany)



- 1 Clinical baseline situation of the vertical and horizontal defect.
- 2 Vertical and horizontal bone defect after surgical opening via a mid-crestal incision.
- 3 3-D printed Titanium scaffold with pre-specified fixation options in the apical area of the "protect" structure.
- 4 Try-in of Yxoss CBR® protect into the defect region. A bone defect of approx. 6 mm is revealed. Geistlich Biomaterials are prepared for placement.
- 5 The scaffold is fixed with a single FYxoss screw in the crestal area. In accordance with the backward planning one implant could be inserted immediately (pos. 14). In addition an external sinus lift was performed.
- 6 Tension-free closure of the double layer mucoperiosteal flap over Yxoss CBR® protect using resorbable deep mattress and single interrupted sutures.
- 7 Clinical image shows an irritation-free healing after the extensive horizontal and vertical augmentation.
- 8 Yxoss CBR® protect successfully prevented hard and soft tissue ingrowth in the apical area. The bone maturation was not affected.
- 9 Thanks to the Yxoss CBR® protect design the scaffold could be removed more easily after a healing period of 6 months.
- 10 Fully regenerated and matured bone with implants.
- 11 After inserting the two additional implants and the healing abutments the flap is closed once again.
- 12 Radiographic image visualizes the implant positions within the regenerated bone.

# Yxoss CBR®



Customized Solutions

Easy ordering at  
[www.reoss.us](http://www.reoss.us)

**Geistlich Pharma North America, Inc.**  
Princeton, NJ 08540  
Customer Care Toll-free: 855-799-5500  
[yxoss@geistlich-na.com](mailto:yxoss@geistlich-na.com)  
[dental.geistlich-na.com](http://dental.geistlich-na.com)

**Manufacturer of Yxoss CBR®**  
ReOss® LLC  
Echterdingenstraße 57  
D-70794 Filderstadt  
[www.reoss.us](http://www.reoss.us)

marketed by  
**Geistlich**  
the regeneration  
experts



Yxoss CBR® classic



Yxoss CBR® protect

**Yxoss CBR®** is an innovative solution for the regeneration of complex alveolar bone defects by using CBCT data in combination with 3-D printing technology. It has simplified the surgical technique for gaining new bone height and width by providing accuracy of fit, volume stability, and predictability.



## Geistlich Bio-Oss®

Stable scaffold for new bone.<sup>1,2,3,4</sup> The slow resorption of Geistlich Bio-Oss® increases the stability of the augmentation material<sup>5</sup> – the best prerequisite for long-term implant survival rates.<sup>6</sup>



## Geistlich Bio-Gide®

Stabilizes the grafted area and protects bone particles from dislocation for optimal bone regeneration.<sup>7</sup> The natural collagen structure allows homogeneous vascularization, supports tissue integration and wound stabilization.<sup>8</sup>

The combination of flexibility, good adhesion, and tear resistance contribute to easy handling, in turn saving time, and simplifying the surgical procedure.<sup>9</sup>

- 1 Orsini G et al., J Biomed Mater Res, B: Appl Biomater 74B, 2005; 448–57.
- 2 Piattelli M et al., Int J Oral Maxillofac Implants 1999; 14: 835–40.
- 3 Sartori S, et al., Clin Implants Res 2003; 14: 369–72.
- 4 Traini T et al., J Periodontol. 2007 May; 78(5): 955–961.
- 5 Orsini G et al., Oral Diseases. 2007; 19: 357–368.
- 6 Jung R et al., Clin Oral Implants Res. 2013 Oct; 24(10): 1065–73.
- 7 Perelman-Karmon M et al., Int J Periodontics Restorative Dent. 2012 Aug; 32(4): 459–65.
- 8 Rothamel D et al., Clin. Oral Implants Res. 2005; 16(3): 369–378.
- 9 Data on File. Geistlich Pharma AG, Wolhusen, Switzerland.



For more information incl. a comprehensive brochure with both Yxoss CBR® versions please visit:  
[www.reoss.eu](http://www.reoss.eu)  
[www.geistlich-pharma.com](http://www.geistlich-pharma.com)

**CAUTION:** Federal law restricts these devices to sale by or on the order of a dentist or physician. For more information on contraindications, precautions, and directions for use, please refer to the Instructions for Use at: [dental.geistlich-na.com/ifu](http://dental.geistlich-na.com/ifu)