

## Micross Curve


Easily reaches challenging areas

**Micross Curve** represents the next evolution in the Micross family, making minimally invasive bone harvesting even easier.



### Minimize discomfort with extra access

The Micross Curve builds on the principle of Micross. Just like Micross, it was developed for minimally invasive techniques. Thanks to the curved design, it is even easier to use in hard-to-reach areas. Collect only the bone that is absolutely needed and reduce invasiveness – your patients will love it.




## Micross Curve

Ergonomic design, simple removal



### New 27° curvature . . .

The Micross Curve ergonomic design allows the user to easily reach retromandibular, palatal and other challenging areas.



### . . . same easy retrieval

You can enjoy the design advantages of Micross Curve without worrying about how to retrieve the bone shavings. The chamber opens the same way as Micross, simply by opening the blue lock and pushing the chamber forward.

Learn more  
about our solutions:



[www.geistlich.us](http://www.geistlich.us)

Ref.-Nr.	Product	Description
Ref. 5542	Micross Curve	1 disposable unit, chamber capacity: 0.25 cc

CAUTION: Federal law restricts these devices to sale by or on the order of a dentist or physician.

For more information on contraindications, precautions, and directions for use, please refer to the Instructions for Use at: [dental.geistlich-na.com/ifu](http://dental.geistlich-na.com/ifu)